



FRAME 708D T30 B80

Plan No. / Plannr.

E000811/C

Scale / Schaal: 1:5	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd
Date / Datum	11/04/2012	29/03/2017	29/03/2017
Init. / Init.	TT	SV	PJ



Material / Materiaal			
No.	Plan	Article / Artikel	Number
1	X000582	1005740	1
2	X002131	1992422 (L)	1
	X002130	1992421 (R)	1
3	X002132	1992423	1
4	X002163	1488100	1
5	X002169	1992791	1
6	X003596	1998528	1
7	XP00094	2990044	1
8	XP00095	2990045	1
9	X003579	1998414	2
10	YP00343	18111.. (L)	1
	YP00346	18112.. (R)	1
11	Ø2x12 DIN1475	1992597	1
12	fork plastic	1991447	1
13	M8x76 DIN7984/8.8	1998993	1
14	M8	1998994	1
15	M12 DIN127A	1001572	1
16	M6 DIN125A	1001564	2
17	M12 DIN985	1001684	1
18	M12	1003294	1
19	M12x100 DIN931/8.8	1009034	1
20	M6 DIN985	1001680	1
21	M6x12	1008492	1
22	M6x45	1991017	1
23	X005513	1992600	1
24	Ø4.8x22 DIN7981	2990070	2
25	X001293	1400504	1
26	X001294	1400604	1
27	X003662	1998693	1
28	X002437	1992699	1
29	M16 DIN128A	1993237	1
30	XC00414	3991378	1
31	Ø12x107	1905243	1
32	X000203	G000791	1

FRAME 708D T30 B80				Plan No. / Plannr.
				E000811/C
Scale / Schaal: 1:5	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	11/04/2012	29/03/2017	29/03/2017	
Init. / Init.	TT	SV	PJ	