



Material / Materiaal			
No. / Nr.	Plan	Article / Artikel	Quantity / Aantal
1	X001398	1009689	1
2	Y013427	1676300	1
3	Y013685	1685600	1
4	Y013683	1685500	1
5	M6x30 ISO7380	1009588	1
6	M6 DIN125A	1001564	3
7	M6x15	1001577	1
8	M6 DIN985	1001680	1
9	M6x16 ISO7380	1008459	4
10	XP00816 (L)	2991000	1
	XP00817 (M)	2991001	1
11	Ø25 M10x1	1999935	1
12	M6x45 ISO7380	1991017	2

A: Mounting holes seat widths 40-45-50  
 B: Mounting holes seat width 35

1902124

E000443/E

DISMOUNTABLE B22 TYPE INOVYS			Plan No. / Plannr.
B22 WEGKLAPBAAR TYPE INOVYS			E000443/E
Scale / Schaal: 2:5	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd
Date / Datum	03/03/2008	01/10/2019	01/10/2019
Init. / Init.	SD	SD	PJ

