


EP00029/F

ELECTRIC BACKREST FRAME SPRINGER KIDS				Plan No.
1/2				EP00029/F
Scale :	1:4	Drafted	Verified	Confirmed
Date	23/03/2011	04/12/2012	04/12/2012	04/12/2012
Init.	AS	BS	TK	



EP00029 /F

Material			
No.	Number	Article	Plan
1	YP02544	2803354	1
2	YP02540	2803254	1
3	YP03008	28222..	1
4	---	1003299	1,2m
5	---	1906790	1
6	YP02160	1893554	2
7	YP02159	1893454	2
8	XP00550	2990635	2
9	M6x12	1007000	4
10	M6x16	1008459	4
11	DIN7991-M6x35	2991297	4
12	DIN125A-M6	1000973	4
13	DIN985-M6	1001680	4
14	DIN7991-M8x45	2991298	2
15	M10x50	1009335	1
16	DIN125A-M10	1001558	2
17	DIN985-M10	1001683	1
18	X000862	1006463	2
19	DIN985-M8	1001682	2
20	M8x15	1006661	1
21	blue	2991367	1
	orange	2991368	1
22	M10x40	1991481	1

ELECTRIC BACKREST FRAME SPRINGER KIDS				Plan No.
2/2				EP00029/F
Scale : 1:4	Drafted	Verified	Confirmed	
Date	23/03/2011	04/12/2012	04/12/2012	
Init.	AS	BS	TK	